



OMI HORSES Pulsed Electromagnetic Therapy



Pulsed Electromagnetic Therapy Equine Product Line

Description

OMI's new Equine Product Line offers a full range of PEMF therapy products that are developed to combat the stresses inflicted on the horse during equestrian sports. The devices work in the frequency range from 1 to 99 Hz and it has three different programmes. Powered by AC power supply. Material: PU Leather.

Recommendation

With OMI PEMF Therapy, your horse can receive the ideal therapy suited for all tissue types including muscles, ligaments and tendons to promote healthy joint function and optimize biomechanical efficiency.

Advantages

All equine products come with adjustable straps to make the applicator fit properly. Simple to use, no training required. Durable and easy to clean. We offer a 2-year warranty for our Equine Product Line.

Instructions for use

The product is designed to be easy to use with no professional experience required. Due to the low intensity of PEMF your horse will not experience any discomfort.

Three different programmes are available: P1 (default programme, the frequency gradually increases from 3 Hz to 38 Hz and then drops back), P2 (default programme, the frequency continually and randomly changes between 1–99 Hz), P3 (the frequency is adjustable between 1 and 99 Hz).

How it works

The gentle magnetic fields produced by our device induce electrical current inside your horse's body that helps to restore the blocked ion flow and stimulates regeneration processes.

Packaging contains

OMI HORSES PEMF Applicator + controller + power supply + user manual

Package Offer: All 5 types of applicators + 5 controllers + 5 power supplies + user manuals

Technical information

Country of Origin: European Community

Delivery Term, Currency and Price: [ask for quote](#)

Payment Terms: on delivery

Pick-up Point: Budapest

Minimum Order Quantity: [ask for quote](#)

Product Range



PEMF WRAP for Rear Leg



PEMF WRAP for Front Leg



PEMF WRAP for Neck



PEMF Blanket



PEMF BAND for Shoulder

product category: wellness product – equine product line

website: www.omihorses.com

www.oxfordmedicals.com/pemf-for-horses.html



oh-pis-gb-180313 • v1.5

Manufactured by

Oxford Medical Instruments Limited

2nd Floor Curzon House, 24 High Street, Banstead, Surrey, England, SM7 2LJ

www.oxfordmedicals.com